

REPOSE CAN BE CLEANED AND REUSED



HAND WASH



**DO NOT
DRY CLEAN**



**DO NOT
TUMBLE DRY**



**DO NOT
IRON**

**REPOSE CAN BE CLEANED WITH SOAP OR DETERGENT
AND WATER. REPOSE CAN BE DISINFECTED IN
ACCORDANCE WITH INFECTION CONTROL POLICY FOR
STANDARD MATTRESSES.**

**FOR FURTHER INSTRUCTIONS, SEE PUMP AND
INSTRUCTION LEAFLET**

© REPOSE WORLD WIDE PATENT PROTECTION

**FLAME RETARDANT - COMPLIES WITH BS7175:1989
IGNITION SOURCE (CRIB 5)**



THIS SIDE UP

